

Logs

The log file section provides access to the raw log files produced by the server. This page always displays log files available in a reverse chronological order so the latest log file is always on top. You can change the number of log files presented using the *Show Top* textbox.

There are two types of log files available:

- PMA.core log: Those are PMA.core's internal logging system with information about *operations on the server, exceptions, errors, access log* etc.
- IIS log: Those are raw IIS log files, containing all the requests made to the server along with all IIS logged info like *time taken, host, http status code* etc

| Available log files | | |
|---|--|--------------------------------|
| Home > Log files | | |
| <div> <div>Top 20 most recent files</div> <div> <div>Show top</div> <div>20</div> <div></div> </div> </div> | | |
| Date | PMA Core log | IIS log |
| Tuesday, March 29, 2022 | tile_server_log-2022-03-29.log | u_ex220329.log |
| | tile_server_log-2022-03-28.log | u_ex220328.log |
| Monday, March 28, 2022 | tile_server_log-2022-03-27.log | u_ex220327.log |
| Sunday, March 27, 2022 | tile_server_log-2022-03-26.log | u_ex220326.log |
| Saturday, March 26, 2022 | tile_server_log-2022-03-25.log | u_ex220325.log |
| Friday, March 25, 2022 | tile_server_log-2022-03-24.log | u_ex220324.log |
| Wednesday, March 23, 2022 | tile_server_log-2022-03-23.log | u_ex220323.log |

More insight into errors is provided via [Elmah](#).

From:
<https://docs.pathomation.com/pma.core/2.0.2/> - **PMA.core 2.x**

Permanent link:
<https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=logs&rev=1649063016>

Last update: **2022/04/04 12:03**



