


Submitted data




The data submitted by each user for each form can be obtained in the “Submitted data” tab in the formats:

- HTML
- CSV
- XML
- ARFF (Attribute-Relational file format)

Forms

Home > Forms > Submitted data

 Submitted data

NSCLC form					
Version	Results				
1	2170				Arff

HTML

CSV

XML

ARFF

From: <https://docs.pathomation.com/pma.core/2.0.2/> - **PMA.core 2.x**

Permanent link:
https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=submitted_data&rev=1648474500

Last update: **2022/03/28 16:35**

