

System requirements

System requirements, software prerequisites, and other recommendations are described in a separate document "[System Requirements](#)" (but probably not as elaborate as in that document)

On-premise installation

Component	Minimum	Suggested
CPU cores	4	16
RAM	8GB	16GB
Disk	250GB HDD	500GB SSD
Cache disk	100GB HDD	1TB HDD
Network interface	100MBps	1Gigabit

* SQLServer on same machine or otherwise? * CPU / Number of cores * Memory * SQLite requirements?

Additional comments regarding SQLServer

SQL CE vs. "full" SQL Server

Installation on an AWS Virtual Machine

* EC2 requirements, * RDS requirements, * store slides on local hard disk (additional partition configured on the EC2 instance) on S3 bucket storage?

Installation on an Azure Virtual Machine

* VM requirements * SQLServer requirements * store slides on local hard disk (additional partition configured on the EC2 instance) on Azure (blob/table/data lake?) storage?

Additional information on cache storage requirements

From:
<https://docs.pathomation.com/pma.core/2.0.2/> - PMA.core 2.x

Permanent link:
https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=system_requirements&rev=1648020172

Last update: 2022/03/23 10:22

