

# User states

A PMA.core user can be in one or more states defined here

## Administrator

The user is an Administrator, which has access to the Administrator UI of PMA.core as well access to all admin API features

## Suspended

The user has suspended all operations, he cannot login or perform any API operation. This is as good as if a user deleted.

## Locked out

The user has been locked out of the system because of [too many incorrect attempts](#) at login or by an administrator

From:

<https://docs.pathomation.com/pma.core/2.0.2/> - **PMA.core 2.x**

Permanent link:

[https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=user\\_states&rev=1648808581](https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=user_states&rev=1648808581)

Last update: **2022/04/01 13:23**

