


Elmah log

Elmah error logs are a standard way of capturing runtime errors in software. While these logs are not specific to root directories, they can help troubleshoot intermittent problems. Examples include:

- connection information is (temporarily) no longer correct
- a file format may have been modified by its vendor and now no longer renders correctly



Host	Code	Type	Error	User	Date	Time
ASTRA4022	0	WebResourceAccessDenied	Log on failure		2/1/2022	6:14 AM
ASTRA4022	0	WebResourceAccessDenied	Log on failure		2/1/2022	6:16 AM
ASTRA4022	0	Http	Threadpool (Process ID 300) ran out of resources with another process and has been chosen as the default victim. Restart the transaction.		2/1/2022	12:28 PM
ASTRA4022	0	WebResourceAccessDenied	Log on failure		2/1/2022	6:14 AM
ASTRA4022	0	Web	The remote server returned an error (403): Forbidden.		2/8/2022	7:48 AM
ASTRA4022	0	Web	The remote server returned an error (403): Forbidden.		2/8/2022	7:48 AM
ASTRA4022	0	Web	The remote server returned an error (403): Forbidden.		2/8/2022	6:38 AM
ASTRA4022	0	Web	The remote server returned an error (403): Forbidden.		2/8/2022	6:38 AM
ASTRA4022	0	Web	The remote server returned an error (400): Bad Request.		2/8/2022	6:38 AM

From:

<https://docs.pathomation.com/pma.core/2.0/> - **PMA.core 2.x**

Permanent link:

<https://docs.pathomation.com/pma.core/2.0/doku.php?id=elmah&rev=1647256401>

Last update: **2022/03/14 14:13**

