

Instances

This section is particular relevant in the context of [geo-replication scenarios](#).

Instances

Home > Instances

Instances Add instance

Name	Description	Current instance	
Brazil		☐	Edit
California		☐	Edit
Frankfurt		☒	Edit
Mumbai		☐	Edit

Showing 1 to 4 of 4 entries

Previous 1 Next

Set current
Delete

Adding a new instance

Add instance

Home > Instances > Create

Add instance

Name

Description

[Create](#) [Back to List](#)

From:
<https://docs.pathomation.com/pma.core/2.0/> - **PMA.core 2.x**

Permanent link:
https://docs.pathomation.com/pma.core/2.0/doku.php?id=settings_instances&rev=1648487052

Last update: **2022/03/28 20:04**

