


## Submitted data

The data submitted by each user for each form can be obtained in the “Submitted data” tab in the formats:

- Unordered List ItemHTML
- Unordered List ItemCSV
- Unordered List ItemXML
- Unordered List ItemARFF (Attribute-Relational file format)



The screenshot shows a web interface titled "Forms". Below the title is a breadcrumb trail: "Home > Forms > Submitted data". There are two tabs: "Submitted data" (active) and another tab with a grid icon. Below the tabs is a table for the "NSCLC form". The table has columns for "Version" and "Results", and then six columns with icons representing different data formats: HTML, CSV, XML, ARFF, and two others. The data row shows Version "1" and Results "2170".

| NSCLC form |         |                                                                                     |                                                                                     |                                                                                       |                                                                                       |      |
|------------|---------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------|
| Version    | Results |                                                                                     |                                                                                     |                                                                                       |                                                                                       |      |
| 1          | 2170    |  |  |  |  | Arff |

From:

<https://docs.pathomation.com/pma.core/2.0/> - **PMA.core 2.x**

Permanent link:

[https://docs.pathomation.com/pma.core/2.0/doku.php?id=submitted\\_data&rev=1591180417](https://docs.pathomation.com/pma.core/2.0/doku.php?id=submitted_data&rev=1591180417)

Last update: **2020/06/03 13:33**

