





## Submitted data

The data submitted by each user for each form can be obtained in the “Submitted data” tab in the formats:

- HTML
- CSV
- XML
- ARFF (Attribute-Relational file format)



The screenshot shows a web interface titled "Forms". Below the title is a breadcrumb trail: "Home > Forms > Submitted data". A tab labeled "Submitted data" is selected. Below this, a section titled "NSCLC form" contains a table with the following data:

Version	Results				
1	2170				Arff

From:

<https://docs.pathomation.com/pma.core/2.0/> - **PMA.core 2.x**

Permanent link:

[https://docs.pathomation.com/pma.core/2.0/doku.php?id=submitted\\_data&rev=1591180478](https://docs.pathomation.com/pma.core/2.0/doku.php?id=submitted_data&rev=1591180478)

Last update: **2020/06/03 13:34**

