




# Submitted data

The data submitted by each user for each form can be obtained in the “Submitted data” tab in the formats:

- [HTML](#)
- [CSV](#)
- [XML](#)
- [ARFF](#) (Attribute-Relational file format)

Forms						
Home > Forms > Submitted data						
Submitted data						
NSCLC form						
Version		Results				
1		2170				 Arff

**HTML**

**CSV**

**XML**

**ARFF**

From:  
<https://docs.pathomation.com/pma.core/2.0/> - **PMA.core 2.x**

Permanent link:  
[https://docs.pathomation.com/pma.core/2.0/doku.php?id=submitted\\_data&rev=1648474500](https://docs.pathomation.com/pma.core/2.0/doku.php?id=submitted_data&rev=1648474500)

Last update: **2022/03/28 16:35**

