

User monitoring

Overview table

User details view

By clicking on the user's login in the overview table, you can see a summary of the user's settings and activity. This page provides an overview of the user's information and the root directories they can access.

The screenshot shows a web application interface for 'User admin Details'. The page has a breadcrumb trail: Home > Users > Details. There are tabs for Summary, Edit, Activity log, Audit trail, and ACL. The 'Summary' tab is active.

User account details

First Name	Administrator	Created by	System
Last Name	Pathomation	Created on (UTC)	4/15/2019 12:00:00 AM
Email	admin1@pathomation.com	Last modified by	
Administrator	<input checked="" type="checkbox"/>	Last modified on (UTC)	
Suspended	<input type="checkbox"/>	Password	Remind password
Session ID	lmkX1mboJl_lo8EAeghZPQ2 User is online		

ACL summary (access control lists)

Public root-directories	_sys_ref	_sys_aes_s3	hgx	hgx_xdr
	invalid_local	invalid_s3_cred	invalid_s3_bucket	Test
Private root-directories	No private root-directories defined.			

Access rights

Activity log

The activity log section provides a detailed list of events that involve the displayed user.

1227 x 817 px (displayed at 772 x 514)

User admin Log

Home > Users > Log

Summary Edit Activity log Audit trail ACL

System Access Log

Session ID	First Accessed	Last Accessed	Host Address	Host Name	Application	Success	Failure Reason
ImkX1mboJl_lo8EAeghZPQ2	4/15/2019, 9:50:49 AM	4/15/2019, 2:15:31 PM	85.73.248.94	85.73.248.94	PMA core	true	
	4/15/2019, 4:34:35 AM	4/15/2019, 4:34:35 AM	127.0.0.1	127.0.0.1	test	false	Invalid username or password
upjGsSqqVKQpkpbK268IA2	4/15/2019, 4:32:52 AM	4/15/2019, 9:45:22 AM	85.73.248.94	85.73.248.94	PMA core	true	

Showing 1 to 3 of 3 entries

Previous1Next

Slide Access Log

First Accessed	Last Accessed	Session ID	Slide Path	Slide Name
4/15/2019, 12:49:17 PM	4/15/2019, 2:15:32 PM	ImkX1mboJl_lo8EAeghZPQ2	._sys_ref/brightfield/Smartzoom	cmu-1.sz
4/15/2019, 12:48:25 PM	4/15/2019, 2:15:32 PM	ImkX1mboJl_lo8EAeghZPQ2	._sys_aws_s3/brightfield/Hamamatsu/CMU-1-VMS	CMU-1-40x - 2010-01-12 13:24:05 vms
4/15/2019, 12:47:59 PM	4/15/2019, 2:15:32 PM	ImkX1mboJl_lo8EAeghZPQ2	._sys_ref/brightfield/Hamamatsu/CMU-1-VMS	CMU-1-40x - 2010-01-12 13:24:05 vms
4/15/2019, 12:47:57 PM	4/15/2019, 2:15:32 PM	ImkX1mboJl_lo8EAeghZPQ2	._sys_aws_s3/brightfield/Generic_Png	png01_julia_4K.png
4/15/2019, 12:47:54 PM	4/15/2019, 2:15:32 PM	ImkX1mboJl_lo8EAeghZPQ2	._sys_ref/brightfield/Generic_Png	png01_julia_4K.png

Audit trail

The audit trail section displays a log of modification events that involve the current user.

1230 x 471 px (displayed at 854 x 327)

Audit trail for user: test

Home > Users > Audit trail

Summary Edit Activity log Audit trail ACL

Audit trail

EventId	EventId	ActionType	EventDate	UserLogin	UserHost	Administrator	CanAnnotate	Description	Email	FirstName	LastName
0	174	Update	4/15/2019, 4:34:15 PM	admin	85.73.248.94	False	True	Test user	test123@pathomation.com	Test	test
0	173	Create	4/15/2019, 4:34:00 PM	admin	85.73.248.94	False	False	Test user	test123@pathomation.com	Test	test

Showing 1 to 2 of 2 entries

Previous1Next

From:
<https://docs.pathomation.com/pma.core/2.0/> - PMA.core 2.x

Permanent link:
https://docs.pathomation.com/pma.core/2.0/doku.php?id=user_monitoring&rev=1646672306

Last update: 2022/03/07 19:58

