

Global settings (administration)


This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

Global settings has 6 fields:

1. Language
2. Application IconNo file chosen
3. Email For Help
4. Conference Server Url
5. Snapshot Filename Template

Annotation settings has 5 fields:

1. Show Annotation Dialog
2. Auto Save Annotations
3. Restart Annotations Drawing
4. Show Slide Info Button
5. Show Full Slide Path In Viewport



From: <https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link: https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=admin_global&rev=1680811111

Last update: **2023/04/06 22:58**

