

























Shared links (administration)

This page list down all the shared links created for a slide or slide folder. It provides detail about each of the shared links that include link type, is password protected and expiry date. You can also delete or suspend a shared link.

Shared Links							
Id	Name	Comment	Type	Is protected by password	Expires At	Actions	
03U6dma			ROI	Not Protected	Jan 01 0001 00:00 Expired		
04qVAKUj	testje Cleg		Slide	Not Protected	Never		
05mKU1UQ			Slide	Not Protected	Never		
0UGVFUJ			Grid	Not Protected	Never		
0JlmmT			Slide	Not Protected	Never		
0SHKCT7aY	Cleg test		ROI	Protected	Jul 08 2021 21:00 Expired		
0wZ6CH9			ROI	Not Protected	Never		
1Aqpsdq			ROI	Not Protected	Never		
113s0HMQ			Slide	Not Protected	Aug 17 2020 19:00 Expired		
1NC933W			Slide	Not Protected	Never		
1KPhwS1E			Grid	Not Protected	Never		
1M3F0My			Slide	Not Protected	Never		

From:

<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:

https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=admin_links

Last update: **2023/04/06 22:57**

