


Shared links (administration)

This page list down all the shared links created for a slide or slide folder. It provides detail about each of the shared links that include link type, is password protected and expiry date. You can also delete or suspend a shared link.



| ID | Name | Summary | Type | Is protected by password | Expires at | Actions |
|------------|------------|---------|------|--------------------------|-------------------|--------------------|
| 1000000001 | 1000000001 | | URL | Not Protected | Jul 31 2024 11:59 | [Suspend] [Delete] |
| 1000000002 | 1000000002 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000003 | 1000000003 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000004 | 1000000004 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000005 | 1000000005 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000006 | 1000000006 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000007 | 1000000007 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000008 | 1000000008 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000009 | 1000000009 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000010 | 1000000010 | | URL | Not Protected | None | [Suspend] [Delete] |
| 1000000011 | 1000000011 | | URL | Not Protected | None | [Suspend] [Delete] |

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=admin_links&rev=1648661905

Last update: **2022/03/30 20:38**

