

Shared links (administration)

This page list down all the shared links created for a slide or slide folder. It provides detail about each of the shared links that include link type, is password protected and expiry date. You can also delete or suspend a shared link.

ID	Name	Summary	Type	Is protected by password	Expires at	Actions
10000001	10000001		URL	Not Protected	Jul 31 2024 11:59	Delete
10000002	10000002	Sample link	URL	Not Protected	None	Delete
10000003	10000003		URL	Not Protected	None	Delete
10000004	10000004		URL	Not Protected	None	Delete
10000005	10000005		URL	Not Protected	None	Delete
10000006	10000006	Sample link	URL	Protected	Jul 31 2024 11:59	Delete
10000007	10000007		URL	Not Protected	None	Delete
10000008	10000008		URL	Not Protected	None	Delete
10000009	10000009		URL	Not Protected	None	Delete
10000010	10000010		URL	Not Protected	None	Delete
10000011	10000011		URL	Not Protected	None	Delete

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=admin_links&rev=1648661905

Last update: **2022/03/30 20:38**

