


Ribbon layout (administration)

This admin page allows you to setup default Ribbon layout that consists of different tabs and controls. Each tab consists of different controls associated to it and you can configure appearance of tabs and controls using this page.



An [exhaustive list of all available commands](#) is available elsewhere on this site.

Below are some of the tabs and corresponding controls used in PMA.Studio.

Home tab has 8 sections:

1. View
2. Grid
3. Snap
4. Magnify
5. Directory
6. Slide
7. Share
8. Search



Grid tab supplies instruments to manage grid and it has 6 sections:

1. ItemSize
2. Synchronize
3. Rows/Columns
4. Slide

- 5. Snapshot
- 6. Share



Annotations tab (available in PMA.Studio) supplies instruments for working with annotations and has 4 sections:

- 1. Annotations
- 2. Shapes
- 3. Export
- 4. Import



Pre-set Annotations tab (available in PMA.Studio) consists of 3 sections:

- 1. Annotations
- 2. Stamps
- 3. Protocols



Conferencing tab has 2 sections:

- 1. Conferencing
- 2. Server



Configure tab has 4 sections:

- 1. Server
- 2. Panels
- 3. Pre-set layouts
- 4. Viewport



From: <https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link: https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=admin_ribbon&rev=1674470758

Last update: **2023/01/23 13:45**

