


Simple shapes

Below are the set of simple shapes that can be added to a annotation.



1. Point
2. Circle
3. Line
4. Arrow
5. Ellipse
6. Rectangle
7. Text
8. Landmark 1
9. Landmark 2
10. Landmark 3

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=annotations_simple&rev=1650554609

Last update: **2022/04/21 18:23**

