


Version retrieval (automation)

For monitoring purposes, it is useful to determine the PMA.studio version used through a back-end API call. All of Pathomation's platform components have this option, actually.

In order to determine the installed version of PMA.studio, you can issue the call:

<https://server/pma.studio/api/json/getversioninfo>

The result is a JSON string formatted as A.B.C.D, with A, B, C being numeric. D can be either numeric or alphanumeric.



Reporting this version number to Pathomation when [asking for support through our helpdesk](#) is always helpful, too.

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=auto_version&rev=1643105135

Last update: **2022/01/25 13:05**

