

## Slide information panel

The Info panel displays available slide information for currently selected slide. Currently selected slide is highlighted with blue dashed border. Columns are pretty self-descriptive and only thing to mention here is that barcode and overview is clickable and by clicking on that images it will be opened in a new browser tab.

From:

<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:

[https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=panel\\_info&rev=1677772733](https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=panel_info&rev=1677772733)

Last update: **2023/03/02 18:58**

