

## Ribbon

Ribbon is a panel that sticks to the top of the screen and consists of tabs and controls. Appearance of tabs and controls are being controlled by ribbon layout, which can be specified in the Admin area. This paragraph will describe behavior of controls, but not their position, since it can be changed by user.

**Home tab** Home tab has 7 sections:

1. View
2. Grid
3. Directory
4. Slide
5. Share
6. Search
7. Conferencing

**Annotations tab (available in PMA.Studio)** Annotations tab supplies instruments for working with annotations and has 3 sections:

1. Annotations
2. Shapes
3. Export

**Grid tab** Grid tab supplies instruments to manage grid and it has 6 sections:

1. ItemSize
2. Synchronize
3. Rows/Columns
4. Slide
5. Snapshot
6. Share

**Configure tab** Configure tab has 2 sections:

1. Server (available in PMA.Studio)
2. Layout

From:

<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:

<https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon&rev=1593505169>

Last update: **2020/06/30 11:19**

