

Ribbon

Ribbon is a panel that sticks to the top of the screen and consists of tabs and controls. Appearance of tabs and controls are being controlled by ribbon layout, which can be specified in the Admin area. This paragraph will describe behavior of controls, but not their position, since it can be changed by user.

Home tab has 7 sections:

1. View
2. Grid
3. Directory
4. Slide
5. Share
6. Search
7. Conferencing

Annotations tab (available in PMA.Studio) supplies instruments for working with annotations and has 3 sections:

1. Annotations
2. Shapes
3. Export

Grid tab supplies instruments to manage grid and it has 6 sections:

1. ItemSize
2. Synchronize
3. Rows/Columns
4. Slide
5. Snapshot
6. Share

Configure tab has 2 sections:

1. Server (available in PMA.Studio)
2. Layout

From:

<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:

<https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon&rev=1593505202>

Last update: **2020/06/30 11:20**

