

Ribbon

Ribbon is a panel that sticks to the top of the screen and consists of tabs and controls. Appearance of tabs and controls are being controlled by ribbon layout, which can be specified in the Admin area.

File Menu tab

This is the starting tab of PMA.studio and provides set of menu items used by the application. It includes:

1. About us
2. Show EULA
3. Manual
4. Keyboard shortcuts
5. File formats

HOME tab

HOME tab consists of different sections that contain set of controls for:

1. Layout
2. Viewport
3. Navigate
4. Grid
5. Snap
6. Magnify
7. Slide
8. Tray
9. Share



GRID tab

GRID tab consists of different sections that contain set of controls:

1. Grid Size
2. Synchronize
3. Rows and Columns (R/C)
4. Slide Zoom in and Zoom out
5. Snapshot
6. Share



ANNOTATIONS tab

ANNOTATIONS tab consists of different sections that contain set of controls for:

1. Annotations
2. Style

- 3. Select Annotation
- 4. Simple Shapes
- 5. Polygon Shapes
- 6. Export
- 7. Import



PRE-SET ANNOTATIONS tab

PRE-SET ANNOTATIONS tab consists of different sections that contain set of controls for:

- 1. Annotations
- 2. Stamps
- 3. Protocols



CONFERRENCING tab

ANNOTATIONS tab consists of different sections that contain set of controls for:

- 1. Conferencing
- 2. Broadcasting



CONFIGURE tab

CONFIGURE tab consists of different sections that contain set of controls for:

- 1. Server
- 2. Panels
- 3. Pre-set layouts
- 4. Viewport



From: <https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link: <https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon&rev=1651177715>

Last update: **2022/04/28 23:28**

