



# Ribbon

We strive to make the user interface (UI) of PMA.studio (and PMA.vue) as recognizable and easy to use as possible.

The top-part of the window is a toolbar implemented as a ribbon. The ribbon is a standard application implement now that is found in both general purpose applications like Microsoft Office, and even in specialized ones like geographic information systems (GIS).



A special mentioning is warranted for the [ribbon's context-menu](#) in the top left-hand corner of the ribbon: this brings up a separate menu with application-wide settings.



From:  
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:  
<https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon&rev=1672140944>

Last update: **2022/12/27 14:35**

