

Ribbon > Configure tab > Panels group > Toggle tray panel

Button identification




Associated command

panel-toggle-tray

Function description

Brings up the [tray panel](#):



Related functionality

[Ribbon > Home tab](#)

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_config_panels_toggle_tray

Last update: **2023/03/03 16:17**

