

Ribbon > Home tab > Layout group > Custom layout

Button identification



Associated command

panels-dialog

Function description

Brings up a dialog that offers the end-user fine-grained control over what [panels](#) are visible:

✚

Items	Show
Slide browser	<input checked="" type="checkbox"/>
Forms Data	<input type="checkbox"/>
Slide Info	<input checked="" type="checkbox"/>
Annotations	<input checked="" type="checkbox"/>
Tray	<input type="checkbox"/>
Conference	<input type="checkbox"/>
Color Adjustment	<input type="checkbox"/>
Shared Links	<input type="checkbox"/>
Search panel	<input type="checkbox"/>
Bookmarks panel	<input type="checkbox"/>
Add grid	<input checked="" type="checkbox"/>
04310688C0010S.scn	<input checked="" type="checkbox"/>
Add IFrame	<input checked="" type="checkbox"/>

Close

Related functionality

[Ribbon > Home tab](#)

From: <https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link: https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_layout_custom

Last update: **2023/03/03 16:25**



