

Ribbon > Home tab > Layout group > Default layout

Button identification



Associated command

Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on the left and a [viewport panel](#) on the right.

Related functionality

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_layout_default&rev=1675860055

Last update: **2023/02/08 15:40**

