

Ribbon > Home tab > Layout group > Default layout


Button identification




Associated command

Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on a left and a [viewport panel](#) on the right.



Related functionality



Last update: 2023/02/08 15:42 ribbon_home_layout_default https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_layout_default&rev=1675860154

From: <https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link: https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_layout_default&rev=1675860154

Last update: **2023/02/08 15:42**

