

# Ribbon > Home tab > Layout group > Default layout

## Button identification




## Associated command

simple-layout

## Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on a left and a [viewport panel](#) on the right.



## Related functionality

[Ribbon > Home tab](#)

Last update: 2023/02/10 16:31 ribbon\_home\_layout\_default [https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon\\_home\\_layout\\_default&rev=1676035914](https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_layout_default&rev=1676035914)

---

From: <https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link: [https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon\\_home\\_layout\\_default&rev=1676035914](https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_layout_default&rev=1676035914)

Last update: **2023/02/10 16:31**

