

Ribbon > Home tab > Viewport group > Toggle zoom controls

Button identification




Associated command


'viewport-toggle-zoom'

Function description

When enabled, the viewport looks like this:



When disabled, the widget is hidden from view:



Related functionality

[Ribbon > Home tab](#)

From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=ribbon_home_viewport_toggle_zoom&rev=1676044134

Last update: **2023/02/10 18:48**

