

Interacting with slides


PMA.studio provides multiple ways to interact with slides that includes:

1. Working with grids
2. Static snapshots
3. Sharing live views
4. Broadcasting and conferencing

UI provides tabs “GRID” and “CONFERENCING” that allows you to interact with slides and perform different actions on it.

(GRID Tab)

Grid tab provides you set of tools that allows you resizing grid size, taking snapshot of grid and live sharing of slides.



From:
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:
<https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=slides&rev=1650991803>

Last update: **2022/04/26 19:50**

