

Context menu for folders

Local site

	Slidename	Slidesize	FileType	Last Modified
↑	..			
📁	Intel		File folder	
📁	local samples		File folder	
📁	PerfLogs		File folder	
📁	Program Files		File folder	
📁	Program Files (x86)		File folder	
📁	scripts		File folder	
📁	Shared folder		File folder	
📁	Slide Samples		File folder	
📁	SvmCache		File folder	

Initiated via mouse right click within the table view borders, it's composed of one menu item :
"Upload folder to remote site".

Upload folder to remote site :

Has been be discussed in details earlier on this tutorial.

Remote site

No context menu options available yet for folders on remote site connections.

From:
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:
https://docs.pathomation.com/PMA.transfer/doku.php?id=context_menu_for_folders

Last update: **2021/07/08 15:51**

