

Context menu for folders

Local site

	Slidename	Slidesize	FileType	Last Modified
	..			
	Intel		File folder	
	local samples		File folder	
	PerfLogs		File folder	
	Program Files		File folder	
	Program Files (x86)		File folder	
	scripts		File folder	
	Shared folder		File folder	
	Slide Samples		File folder	
	SvmCache		File folder	

Upload folder to remote site

Initiated via mouse right click within the table view borders, it's composed of one menu item : "Upload folder to remote site".

Upload folder to remote site :

Has been be discussed in details earlier on this tutorial.

Remote site

No context menu options available yet for folders on remote site connections.

From: <https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link: https://docs.pathomation.com/pma.transfer/doku.php?id=context_menu_for_folders

Last update: **2021/07/08 15:51**

