


Create new folder (remote)




Click on this button pops up the following dialog window :




Name of the folder to be created should be added on the text field then a click on button “OK” for confirmation. If operation is successful, the folder is added and the table view is automatically refreshed to show the newly created folder, otherwise one of the following windows pops up :

- If the operation is prohibited :



- If operation fails (duplicate folder name, access rights..) :



Last update: 2021/07/08
15:50

create_new_folder_remote https://docs.pathomation.com/PMA.transfer/doku.php?id=create_new_folder_remote

From:
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:
https://docs.pathomation.com/PMA.transfer/doku.php?id=create_new_folder_remote

Last update: **2021/07/08 15:50**

