



Click on button “Remote Directory” opens a tree view window to select folder to navigate to. Action takes places by a click on button “OK”, while “cancel” simply cancels operation. It's worth also noting that the folder blue icon refers to a root directory, while the yellow icon refers to a non root one.

From: <https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link: [https://docs.pathomation.com/pma.transfer/doku.php?id=folders\\_navigation\\_tree\\_remote&rev=1589595289](https://docs.pathomation.com/pma.transfer/doku.php?id=folders_navigation_tree_remote&rev=1589595289)

Last update: **2020/05/16 05:14**

