



## Windows

- Extract file PMA.transfer.zip into a destination a folder
- Navigate to destination folder

Name	Date modified	Type	Size
 PMA.transfer-2.0.3.72.exe	1/24/2020 11:08 AM	Application	9,483 KB
 PMA.transfer.ini	3/17/2019 10:11 PM	Configuration settings	1 KB

- Double click on file PMA.transfer-2.0.3.x.exe

The Zip file is optionally shipped with a PMA.transfer.ini file containing a preconfigured connection data, it's loaded automatically when the app is launched.

## MacOS

From:  
<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:  
<https://docs.pathomation.com/pma.transfer/doku.php?id=installation&rev=1598914769>

Last update: **2020/09/01 01:59**

