


# My sites

The node “My Sites” is the root of the site manager tree, it cannot be deleted or renamed.



- Click on button “New Site” creates a new site inside “My Sites”.
- Click on button “New Folder” creates a new folder inside “My Sites”.

From:

<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/pma.transfer/doku.php?id=my\\_sites](https://docs.pathomation.com/pma.transfer/doku.php?id=my_sites)

Last update: **2021/07/08 15:44**

