


- Double click to navigate into a folder

- Double click on this row  (to the top) to navigate to a parent folder

- | | Slidename ▲ | Slidesize | FileType | Last Modified |
|--|-------------|-----------|----------|---------------|
|--|-------------|-----------|----------|---------------|

 Click on column headers to sort content in the table view by slides' name, slides' size, file type (folder or slide) and last modification date in ascendant or descendant order.

These options are exactly the same whether it's a local or remote site connection.

From:
<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:
<https://docs.pathomation.com/pma.transfer/doku.php?id=navigation&rev=1599278482>

Last update: **2020/09/05 07:01**

