

### New features :

- Bypass proxy
- Delete directories on server side (only if empty and not a root one)
- Rename directories on server side (only if empty and not a root one)
- Delete slides on server side
- Rename slides on server side
- Add copy full path option both locally and remotely

### General improvements :

- Directly download to selected folder on local site
- Keyboard input interaction
- New icon for root directories
- Better transfer speed metrics
- Better error messages
- Enhanced integrity check
- Adding counters to all bottom bar buttons

### Bug fixes :

- Refresh table view for local site following successful downloads

From:

<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/pma.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.0.3&rev=1589520834](https://docs.pathomation.com/pma.transfer/doku.php?id=what_s_new_in_version_2.0.3&rev=1589520834)

Last update: **2020/05/16 05:14**

