

### New features :

- User quota enforcement

### General improvements :

- Error handling for uploads to Azure directories
- Improved error handling for server to server transfers
- Support for My Pathomation “add slide” API endpoint & “caller ID” parameter for slides uploaded via PMA.transfer
- Support for Keyboard key <Enter>

### Bug fixes :

- Fixes on the server to server progress bar

From:

<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/pma.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.2.1&rev=1612319011](https://docs.pathomation.com/pma.transfer/doku.php?id=what_s_new_in_version_2.2.1&rev=1612319011)

Last update: **2021/02/03 05:23**

