

What's new in version 2.2.3?

PMA.transfer 2.2.3.9d6e296f92 is an improved version of the previous PMA.transfer with fixed bugs and new features in functionality.

General improvements

- Upload slides and folders to the Amazon s3 cloud.
- Upload slides and folders recursively.
- Displaying information about slides in the cloud directories window.
- Adding a new folder in the cloud directory.

New features in functionality

- Added the possibility to fill folders with slides and subfolders recursively, this means that you can load a folder containing subfolders containing both folders and slides.
- Added the possibility to unload multiple selected slides.
- Added the possibility to delete several selected slides from s3 the cloud directory at once.
- Added the possibility to create a new folder in the cloud S3 directory.

Bug fixed :

- Bug fixed preventing more than one slide from being uploaded.
- Bug fixed preventing making a new folder in the cloud directory.

From:
<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:
https://docs.pathomation.com/pma.transfer/doku.php?id=what_s_new_in_version_2.2.3&rev=1662629474

Last update: **2022/09/08 12:31**

