

# What's new in version 2.2.4.d67f7abd23?

PMA.transfer 2.2.4.d67f7abd23 is an improved version of the previous PMA.transfer with fixed bugs and new features in functionality.

## General improvements

- In this build version upload and download slides and folders execute whru Maven dependency(2.0.0.117)
- The filling of progressive bars executes thru the callback of the HTTP interceptor.
- Upload and download methods allow executing transfer slides to different storages: clouds(Azure and Amazon) and remote servers (file storage)

## New features in functionality

- Added the possibility to rename slide on local PC
- Added the possibility to rename folder on local PC

## Bug fixed

- Fixed progressive bar filling calculation system for slides with the multipart format as .vsi and .mrxs

## Slide formats

**Slide formats with PMA.transfer has been tested:**

- .mrxs
- .svs
- .ndpi
- .vsi
- .jpg
- .png
- .zif
- .tiff

From:

<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.2.4](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.4)

Last update: **2023/02/27 16:06**

